

THE LATEST IN INFUSION TECHNOLOGY

# INTRODUCING THE "Z-PUMP"

Z-800 INFUSION PUMP



Ask about our new 45° mounting bracket, alarm light extension and protocol features.

**THE Z-800 INFUSION PUMP FEATURES:**

- A SIMPLE, EASY TO USE DESIGN
- DURABLE / RELIABLE CONSTRUCTION
- ENHANCED PATIENT SAFETY
- LOWER COST OF OWNERSHIP
- UPGRADEABLE SOFTWARE

Available Through  
**OMNI**  
Medical Supply, Inc.

1-800-860-6664

Designed specifically for nurses in alternate site environments where the pump's simplicity and low cost of operation are ideally suited to basic infusion work flow needs.

**Simple programming and pump operation require less time for training, more time for patient care**

- Simple, intuitive user interface (simplicity=safety)
- Requires the fewest programming steps of any similar pump on the market
- Large, easy-to-read display provides all pump information at a glance

**Built-in safety features protect patients**

- Pump-based free-flow protection; additional set-based free-flow protection also available
- Tamper resistant: Quick one-touch lockout feature
- Alarms include occlusion and air-in-line

**Designed with the bottom line in mind**

- The Z-800 offers the best value infusion pump on the market today with lowest total cost of ownership
- Lower maintenance costs due to durable design
- Additional long-term cost savings achieved over life of the pump with Zyno administration sets

**Ergonomic design is easy to transport, built to last**

- Compact size, light weight, cast aluminum casing and robust design features mean the Z-800 is built to last, even in the most punishing environments
- Pump durability translates into significantly lower maintenance costs
- Compact size and light weight means it is easy to transport for both nurses and patients
- Quiet operation

**Flexible operation for all your needs**

- Multi-Therapy Modes:
  - Continuous Rate/Volume — TPN
  - Continuous Time/Volume — Protocol
  - Ten Step Sequencing — Intermittent

**Upgradeable architecture platform designed for easy future feature enhancements.**

# Z-800



# ZYNO MEDICAL Z-800



## Specifications



**Toll Free: 1-800-860-6664**  
**Fax: 1-248-360-9375**  
**Contact your OMNI**  
**Sales Representative**  
**for more**  
**Information on our**  
**Convenient Programs**  
**For**  
**\*Purchase**  
**\*Rental**  
**\*Leasing**  
**and**  
**\*Usage Agreements**

<b>Mechanism</b>	Linear Peristaltic
<b>Flow Rate Range</b>	1-999 mL/hr
<b>Volume To Be Infused (VTBI)</b>	1-9,999 mL
<b>Accuracy</b>	+/- 5%
<b>Infusion Modes</b>	Continuous Rate/Volume, Continuous Time/Volume, Ten Step Sequencing, TPN, Protocol, Intermittent
<b>Free Flow Protection</b>	Pump, and Optional Set-Based Protection
<b>Housing</b>	Cast Aluminum
<b>Weight</b>	5.7 Pounds (2.6 kg)
<b>Size</b>	8.6 x 5.7 x 5.3 in (H/W/D)
<b>Power</b>	110-240 VAC, 50-60Hz
<b>Leakage Current</b>	<100uA
<b>Battery</b>	9.6V, 4.7Ah NiMH (rechargeable)
<b>Battery Operating Time</b>	8 Hours at 125mL/hr
<b>Air In-Line</b>	Ultrasonic. 14 Level Adjustable
<b>Occlusion Pressure</b>	4-30 PSI Adjustable
<b>Tubing</b>	Zyno Needleless, Non-DEHP Standard Tubing
<b>Alarms</b>	Air-in-Line, Occlusion, Door Open, Battery Empty
<b>Warnings</b>	Infusion Complete - KVO, Pump Unattended, Near End, Low Battery



Visit our web site at [www.omnimedicalsupply.com](http://www.omnimedicalsupply.com)

# Z-800